

TECHNICAL SPECIFICATIONS

Powder free disposable nitrile exam gloves



*

* Image for representational purpose only

Designation	Powder free disposable nitrile exam gloves
Origin	China
Manufacturer	Titanfine®, a brand of Titans Group
SAP Code	HYG_053_S / HYG_053_M / HYG_053_L / HYG_053_XL
EAN Code	NC
Material	Powder-free nitrile
Available sizes	S, M, L, XL
Use	Protection against biological and virucidal agents, accidental contact with chemicals or substances, chemical resistance with concentrated detergent, detergent, hydrocarbon, ether and freon. It is suitable for examination, therapy, diagnostic, laboratory, cleaning, chemical industry, workshops, electronics and food processing use. It is suitable for hospital and ambulatory use.
Standards	Medical Devices Regulation, MDR, Regulation (EU) 2017/745 P.P.E.Cat.III in compliance with Dir. (EU) 2016/425 Hole absence for AQL 1.5 in conformity with UNI EN 455, Part I Physical properties in conformity with UNI EN 455, Part II Biological evaluation in conformity with UNI EN 455, Part III Conformity with UNI EN 374-1-2-4 Conformity with UNI EN 420 - UNI EN 388 Conformity with UNI EN 374-5 (virucidal) Products produced in factories with EN ISO 9001 and EN ISO 13485 certifications.
Expiry date	3 years (from the manufacture date).
Storage conditions	Protect from freezing and avoid excessive heat. Keep dry. Products should be protected from direct sunlight, fluorescent lighting, x-rays, moisture and ozone.
Packaging	100 gloves per box (or 50 pairs) - 10 boxes per carton Box dimensions : 235mm x 125mm x 70m Carton dimensions : 250m x 260mm x 360m

Sizes	S	M	L	XL
Thickness (mm)	0,08 mm	0,08 mm	0,08 mm	0,08 mm
Width (mm)	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	≥ 110 mm
Length (mm)	240 mm	240 mm	240 mm	240 mm
Weight (gr)	3,5 ± 0,3 mm	4,0 ± 0,3 mm	4,5 ± 0,3 mm	5,0 ± 0,3 mm

Additional documents available on request.